

www.sumerinc.com



About Us

Sümer A.Ş. was established in 1981 in Ankara to provide services in the medical device sector. It has aimed advancement since the day of its establishment by also taking growth and compliance with the contemporary technologies and protecting the environmental conditions. Our manufacturing takes place in an area of 25.000 m² in Ankara Organized Industrial Zone.

Our Firm is following the innovations in its sector and in abroad through its research and development unit and its application staff with a strong infrastructure of engineers, and is continuing to produce devices it had developed in computer environment based on such innovations with high technology and to contribute their development so as to be most beneficial for the Turkish medicine.

With this purpose, our Firm is strictly **following the “Quality Management” principles and rules from** design of the products to the after-sale servicing.

Our Firm has been currently certificated for compliance with ISO 9001 quality management system, ISO 13485 medical device quality management system certificate and ISO 14001 environment management system certificate and with product certificates under MDD 93/42/EEC Medical Devices Directive and PED 2014/68/EUPressurized Equipment certificate. Furthermore, our steam sterilizers, and disinfection and washing devices have been certificated by the German accredited body.

Our Firm is successfully implementing several projects supported by Tübitak, KOSGEB, and Ankara Development Agency. SÜMER A.Ş. has been awarded with the following:

Ankara Industry Chamber Award for 2012, 2011 R&D Award, 2012 Success Award.

Our Firm possesses the following certificates;

CE Certificates under the following directives:

ISO 9001,

ISO 13485,

ISO 14001 Quality Management System,

MDD 93/42/EEC Medical Devices,

PED 2014/68/EUPressurized Equipment



Our Vision

To make the Sümer brand a global brand to make our Firm remembered first in the sector.

Our Mission

Our main task is to create designs with competition power in the global sense by taking the priorities of the sector into consideration and being respectful to the environment and people and giving the first priority to the wishes and expectations of customers, and also to produce innovative technological medical products by meeting all the national and international legal requirements.

Our Basic Values

- We are bound up with the Medical Ethical Rules,
- We are people-oriented,
- We respect environment,
- We are creative,
- We are customer-oriented,
- We are innovative,
- We are pro-active,
- We believe in the team spirit.

OPERATION TABLE ACCESSORIES

Our focus is pure – continuous innovation, outstanding quality and usability in all our products

SÜMER offers many kinds of accessories to our valued customers for use in different surgeries. All of our products are designed and manufactured to meet the requirements of the relevant ISO and EN standards.



Instrument Table

Ref. No. 8698918585942



Anesthesia Frame

Ref. No. 8698918583580

Specifications for Instrument Table

Length	50cm
Height	47cm
Width	40cm
Weight	2,9kg
Material	Stainless steel

Specifications for Anesthesia Frame

Length	600mm
Height	min 530mm max 780mm
Width	55mm
Weight	2,1kg
Fastening rail (EUR)	10-310
Material	Stainless steel, Aluminum, Brass






Leg Support

Ref. No. 8698918583634

Specifications for Leg Support

Length	350mm
Height	760mm
Width	195mm
Weight	4kg
Material	Stainless steel, Steel, Polyurethane Foam




Head Support

Ref. No. 8698918585966

Specifications for Head Support

Length	59cm
Height	8,5cm
Width	32,5cm
Weight	4,5kg
Fastening rail (EUR)	10-234,10-234-T
Material	Stainless steel,Aluminum, Antistatic fabrik, Foam, Slow recovery



Kidney Elevator

Ref. No. 8698918585973

Specifications for Kidney Elevator

INTERNAL KIDNEY ELEVATOR

$h_{\max} = 12 \text{ cm}$



Waste Bin

Ref. No. 8698918585959

Specifications for Waste Bin

Length	320 mm
Height	530 mm
Width	150 mm
Weight	6 kg
Material	Stainless steel







Big Clamp

Ref. No. 8698918583665

Specifications for Big Clamp

Length	125mm
Height	160mm
Width	65mm
Weight	1,1kg
Fastening rail (EUR)	10-307
Material	Stainless steel, Aluminum, Brass



Cranial Stabilization(Skull Clamp)

Ref. No. 8698918585980

Specifications for Cranial Stabilization

It is used in brain surgeries

Material	Stainless steel, Aluminum, Brass
----------	----------------------------------



Hand Surgery Set

Ref. No. 8698918583603

Specifications for Hand Surgery Set

Length	62,5cm
Height	70-115cm
Width	34cm
Weight	6,7kg
Fastening rail (EUR)	10-390-K
Material	Stainless steel, Aluminum, Backelit, Glass fibre



Arm Support

Ref. No. 8698918583610

Specifications for Arm Support

Length	605mm
Height	140mm (bar for height adjustment)
Width	18cm
Weight	4,3kg
Fastening rail (EUR)	10-380
Material	Stainless steel, Aluminum



Universal Clamp

Ref. No. 8698918583672

Specifications for Universal Clamp

Length	90mm
Height	60mm
Width	40mm
Weight	0,3kg
Fastening rail (EUR)	10-301
Material	Aluminium anodized




Lateral Support

Ref. No. 8698918583641

Specifications for Lateral Support

Length	27cm
Height	53cm
Width	20,5cm
Weight	2,5kg
Material	Stainless steel, Steel, Polyurethane Foam



Shoulder Support

Ref. No. 8698918583627

Specifications for Shoulder Support

Length	27cm
Height	53cm
Width	20,5cm
Weight	2,5kg
Material	Stainless steel, Steel, Polyurethane Foam



Serum Stand

Ref. No. 8698918583573

Specifications for Serum Stand

Height	94,5-149cm
Width	6,5cm
Weight	2,4kg
Fastening rail (EUR)	10-402
Material	Stainless steel, Plastic




Tibia Support

Ref. No. 8698918585997

Specifications for Tibia Support

Height	82cm
Width	30cm
Weight	2,45kg
Material	Stainless steel, Steel, Polyurethane Foam



Traction Set

Ref. No. 8698918583689

Specifications for Traction Set

Orthopedic traction set is suitable for use with all SÜMER operating tables in orthopedic surgeries.

- Made of 304 quality stainless steel
- Leg strainer function by mechanical lever
- Traction boots with up & down movement and rotation capability around its own axis



Body Belt

Ref. No. 8698918583597

Specifications for Body Belt

Length	241cm
Width	12,5cm
Weight	0,3kg
Material	Terylene




Side Safety Bar

Ref. No. 8698918586000

Specifications for Side Safety Bar

Length	630mm
Height	370mm
Width	40mm
Weight	2,8kg
Material	Stainless steel, Aluminum



Widening Section

Ref. No. 8698918585935

Specifications for Widening Section

Length	300 mm
Height	85 mm
Width	125 mm
Weight	2.9 kg
Material	Stainless steel, Aluminium



Small Clamp

Ref. No. 8698918583658

Specifications for Small Clamp

Length	90mm
Height	60mm
Width	40mm
Weight	0,3kg
Fastening rail (EUR)	10-301
Material	Aluminium anodized



Solution Partners

SIEMENS

Schneider
Electric

seko

HANING

OSV

SKF

DIGITAL

WURTH

SICK



www.sumerinc.com

 **sumer**

Address : Başkent OSB Başkent Bulvarı No:81
Malıköy-Sincan/Ankara/Türkiye
Phone : +90 312 394 61 73 +90 312 503 08 21 (pbx)
Fax : +90 312 394 61 95
E-mail : sumer@sumeras.com

 **SZUTEST**
NOTIFIED BODY 2195

kiwa
agency for progress

 CE

 ISO 13485

 ISO 9001

 ISO 14001